

"Erma" Make Hand Refractometer

Hand Refractometer is kind of precision optical instrument.

Through to provide with a direct reading of the scales we can easily know the concentration or refractive index of liquids

- Easy and very convenient to use.
- High quality, accurate testing result

Range : 0-32% Brix
 : 28-62% Brix
 : 45-82% Brix
 : 58-90% Brix
 : 58-92% Brix



Ref.: 624

Abbe Refractometer

A simplified fast operating instrument for measuring refractive index and sugar percentage. The prism assembly is mounted in horizontal position.

For measuring refractive index and suger percentage directly.

Range : Refractive Index 1.3 to 1.7 nD
 Sugar Percent 0 to 95%

Accuracy : Refractive Index 0.001 direct & 0.0001 estimate.
 Sugar percentage 10 on scale and 0.10 by estimate.

Ref.: 622



Magnifier

- Complete with battery.
- Magnification 8x / 30x with artificial illumination.
- Magnifier without illumination.



Ref.: 177