



## Checker Board

Checker board is used to analyse the opacity, hiding power, spreading area/rate etc of the paint film. It is made up of high quality material which is reliable in performance. Checker board cards are also provided in two types which is given below in the table.

Chart No./ Ref No.	Size	Chart Type
201 - C	10" x 6"	Glossy
301 - C	10" x 6"	Matt
201 - F	10" x 8"	Glossy
301 - F	10" x 8"	Matt



## Plain White Charts

Plain white charts are used to produce true colour samples without any interference of black substrate. We provide plain white charts in three types which are given below in the table .

Chart No./ Ref No.	Size	Chart Type
210	10" x 6"	Glossy
211	10" x 4"	Glossy
310	10" x 6"	Matt



## Checker Stickers

Checker stickers are used to conduct and analyse the test of opacity, hiding power, spreading area/rate of the prepared paint sample film and to provide their customer for various reference purposes.

Chart No./ Ref No.	Size	Chart Type
212	25 mm x 25 mm	Glossy

