

## Colour Reader (colour Muse)

Colour Muse works with an app on your mobile phone to scan and match colours. All you need to do is place the colour muse on the surface and press the scan button on the app. The app is free and will show a list of the best colour matches. It can match any colour oriented product. It is very useful for designers, colour professionals, and technicians due to its portability and accuracy for providing a seamless colour communication tool in colour sensing technology. The Colour Muse is designed to fit in your pocket or bag.

**Standard:** ASTM D 2244, ISO 7724

### Raj Make

#### Features:

- Scan a colour for LAB, LCH, RGB, HSB, CMYK values.
- Compare colour & shades with  $\Delta E$  value.
- View & Export facility via text message, email or media.
- Battery Capacity 280 mAh (5,000 Scans) per charge.
- Compatible with Android or IOS operating systems.
- Highly accurate and inbuilt bluetooth module.
- Auto calibration when turned on
- Create and save your own colour palettes.
- Repeatability 0.10  $\Delta E$
- Measurement Size 4mm

#### Package Includes:

Colour Muse, USB cable, manual.

#### Ordering Informations:

Ref No. 780 Colour Reader

HSN Code: 90278090



## Colour Comparator

Colour Comparator is used to match the colour of the liquid samples in a day light illuminating unit against glass standard built in the test disc. The discs are known as Gardner discs having 18 colours with number (1-18). It is widely used in the oil & solvent industry, paint industry, cosmetic pharmaceuticals and chemical industry. It is very easy to use for visual observation. Just rotate the discs for matching the colour of the liquid sample.

It has two nos. of discs

Gardner Discs A Range 1-9

Gardner Discs B Range 10-18

ASTM Discs A Range 0.5-4.0

ASTM Discs B Range 4.5-8.0

Tubes are used to fill the samples for colour matching.

**Standard:** ASTM D 1544, ASTM D 6166, ASTM D 1500, ISO 2049

#### Package Includes:

Colour comparator apparatus, Discs, gardner empty tubes and manual.

#### Ordering Informations:

Ref No. 779/1 Colour Comparator with Gardner Disc A, B

Ref No. 779/2 Colour Comparator with ASTM Disc A, B

Ref No. 779/3 Lovibond (U.K make) Comparator with Gardner Disc A, B or ASTM Disc A, B

HSN Code: 90278090

