

## Pen Type PH Meter

Imported make

### Specification :

- Range : 0 to 14
- Resolution : 0.1 pH
- Accuracy :  $\pm 0.2$  pH
- Calibration : Digital Automatic Calibration ( One Point )
- Electrode : Glass Sensor and reference tube Electrodes
- Housing : Water resistant

Auto-off function, data - hold function and low - battery indicator.

Ref.: 170 / 1



## Digital pH Meter

Model : pH Cal

### Specification :

- Range : 0 to 14 pH
- Resolution : 0.01 pH
- Accuracy : 0.01 pH  $\pm 1$  digit
- ATC : Manual

Supplied with complete Standard Accessories.  
pH electrode, PH4, PH7 Standard solution (Buffer),  
Electrode stand and manual.

Ref : 170



## Conductivity Meter

Imported make

### Specification :

- EC Range : 0-9999  $\mu$ S
- Temp. Range : 0 - 80°C, 32 - 176° F
- Resolution : 1  $\mu$ S
- Accuracy :  $\pm 0.2$  %
- Calibration : Digital Calibration by push button

### Feature

Measures electrical conductivity (EC) and Temperature  
Automatic Temperature Compensation ( ATC )  
Water - resistant housing  
Digital Calibration ( Push button )  
Auto-off function, data - hold function and low battery indicator.

Ref. 171

