

Heating Mantle

Heating Mantle is an instrument used to apply heat to containers. Glassware containers placed in direct contact with heating mantle for heating without increasing the risk of glassware shattering. Heating mantle with the regulator is made up of mild steel with powder coating for durability. It is used for heating organic liquid placed in glass containers like round bottom flask reaction vessels. We provide heating mantle from 100 ml to 50 ltr with a regulator. Temperature range of the heating mantle is ambient to 450°C. It is widely used in the research and development industry, coating industry, oil and wax industry, pharmaceuticals etc.

Raj Make

Volume Capacity: 100 ml, 250 ml, 500 ml
1, 2, 5, 10, 20, 50 litres

Ordering Informations:

Ref No. 656 Heating Mantle

HSN Code: 84799090

