

Adjustable Micrometer Film Applicator

Adjustable Micrometer Film Applicator is used to draw the film according to requirement by adjusting the micrometer head. It is made up of an aluminium frame and blade. The stainless steel blades are easy to install, calibrate and clean. It enhances the accuracy in the thickness of the film. The blade of this applicator can be adjusted with the help of two micrometer heads from 0 - 10,000 microns with an increment of 10 µm and its least count is 10 µm. The drawn down film width is upto 100 mm.

This instrument is suitable for highly viscous solution and thick fluid on a flat surface. It is made up of high grade non corrosive anodised aluminium. These applicators are designed to be convenient for handling, easy to use / clean and ideal for the paint and coating industry, cosmetic industry and steel industry.

Standard: ASTM D 823

Raj Make

Technical Specifications:

Ref No.	Film Width (mm)	Range (µm)	Least count & Increment
150	100 mm	0 - 10 mm	10 µm

Package Includes:

Adjustable micrometer film applicator and manual

Ordering Informations:

Ref No. 150 Adjustable micrometer film applicator

HSN Code: 90318000

